Health Products and Food Branch

Ottawa

Enumeration of generic *E. coli* and other coliform bacteria from all foods using RAPID'*E. coli* 2 chromogenic media

Microbiological Methods Committee
Microbiology Evaluation Division
Bureau of Microbial Hazards, Food Directorate,
Health Products and Food Branch, Health Canada
Postal Locator: 2204E
Ottawa, Ontario K1A 0K9

Contact the Microbiological Methods Committee: hc.mmc-cmm.sc@canada.ca

1. Application

This method is applicable to the rapid detection of generic *E. coli* and other coliform bacteria to determine compliance with the requirements of Sections 4 and 7 of the *Food and Drugs Act* and/or other relevant federal regulations. This method has been validated for use on all foods for enumeration of coliforms and *E. coli*.

2. Description

Coliforms are a non-taxonomic group of bacteria that are exclusively Gram-negative, non-spore-forming facultative anaerobic, rod-shaped bacteria which ferment lactose to produce acid and gas. The group is composed of enteric bacteria such as *Escherichia*, *Citrobacter*, *Klebsiella* and *Enterobacter*. Enumeration of coliform bacteria may be useful to confirm sanitation of water or food processing environments. *E. coli* is found in the intestine of animals and is abundant in human and animal feces. The bacteria can also be isolated from environmental samples. Fecal coliforms are often used as indicators of contamination and the possible presence of pathogenic organisms. The fecal coliform group consists primarily of *E. coli* and some thermotolerant *Klebsiella*. Although these bacteria are generally non-pathogenic, they can cause infection in immunocompromised hosts. Most probable number (MPN) methods for enumeration of coliform bacteria can be laborious and costly. The use of chromogenic substrates in media has led to development of faster and easier methods for detection, differentiation, and enumeration of target bacteria.

3. Principle

RAPID'*E. coli* 2 medium relies on simultaneous detection of 2 enzymatic activities, β -D-glucuronidase (GLUC) and β -D-galactosidase (GAL). The medium contains 2 chromogenic substrates. One substrate is specific to GAL and results in blue coloration of colonies positive for this enzyme and one substrate is specific to GLUC and results in violet coloration of colonies

positive for this enzyme. Coliforms, other than *E. coli* (GAL+/GLUC–), form blue to green colonies while, specifically, *E. coli* (GAL+/GLUC+) form violet to pink colonies. A count of total coliforms can be obtained by adding the number of blue-green colonies and the number of violet to pink colonies.

4. Definition of Terms

See [Appendix A of Volume 1](#).

5. Collection of Samples

See [Appendix B of Volume 1](#).

6. Materials and Special Equipment

Note: The laboratory supervisor must ensure that completion of the analysis described in this method is done in accordance with the International Standards reference “ISO/IEC 17025:2005 (or latest version): General Requirements for the Competence of Testing and Calibration Laboratories”.

Note: It is the responsibility of the laboratory to ensure equivalency if any variations of the media formulations listed here are used (either product that is commercially available or made from scratch). Please forward equivalency data to the [Editor of Compendium of Analytical Methods](#) for consideration of modification of this method.

1) RAPID'*E. coli* 2 agar

- 1.1 Ready-to-use 100 mL x 6 bottles (Bio-Rad catalog # 3555299)
- 1.2 Ready-to-use 200 mL x 6 bottles (Bio-Rad catalog # 3555297)
- 1.3 Dehydrated medium 500g (Bio-Rad catalog # 3564024)

2) Supplies and Reagents

- 2.1 Butterfield's phosphate buffer
- 2.2 Stomacher bag or stomacher bag with incorporated filter (recommended for meat and highly particulate samples)
- 2.3 Sterile pipettes
- 2.4 Sterile petri dishes
- 2.5 Brilliant Green Lactose Bile broth (BGLB)
- 2.6 Levine's Eosin Methylene Blue agar (L-EMB)
- 2.7 Lauryl Sulfate Tryptose (LST) broth
- 2.8 (EC) *E. coli* broth

3) Equipment

- 3.1 Balance - for weighing food samples
- 3.2 Stomacher, blender, vortex or equivalent
- 3.3 Incubator capable of maintaining temperatures of $37 \pm 1^\circ\text{C}$ and $44 \pm 0.5^\circ\text{C}$
- 3.4 Water bath capable of maintaining temperatures of $44 - 47^\circ\text{C}$

Note: It is the responsibility of each laboratory to ensure that the temperatures of the incubators or water baths are maintained at the recommended temperatures. Where 35°C is recommended in text of the method the incubator may be at 35 ± 1.0°C. Similarly, lower temperatures of 30 or 25 may be ± 1.0°C. However, where higher temperatures are recommended, such as 43 or 45.5°C, it is imperative that the incubators or water baths be maintained within 0.5°C due to potential lethality of the higher temperatures on the microorganism(s) being isolated.

7. Procedure

Each sample unit shall be analyzed individually. The test shall be carried out in accordance with the following instructions:

7.1 Handling of Sample Units

7.1.1 In the laboratory prior to analysis, except for shelf-stable foods, keep sample units refrigerated or frozen, depending on the nature of the product. Thaw frozen samples in a refrigerator, or under time and temperature conditions which prevent microbial growth or death.

7.1.2 Analyze sample units as soon as possible after their receipt in the laboratory.

7.2 Preparation for Analysis

7.2.1 Have sterile Butterfield's phosphate buffer diluent prepared.

7.2.2 Clean the surface of the working area with a suitable disinfectant.

7.3 Preparation of sample

To ensure a representative analytical unit, agitate liquids or free flowing materials until the contents are homogeneous. If the sample unit is a solid, obtain the analytical unit by taking a portion from several locations within the sample unit.

7.3.1 Prepare a 1 in 10 dilution of the food by adding aseptically 11 (50) g or mL to 99 (450) g or mL of sterile Butterfield's phosphate buffer diluent.

Note: Weight or volume in brackets indicates alternate procedure for making dilutions.

7.3.2 Blend, stomach or vortex as required for thorough mixing.

7.3.3 Prepare serial dilutions as necessary in sterile Butterfield's phosphate buffer diluent.

7.4 Enumeration of coliforms and *Escherichia coli*

7.4.1 Plating

7.4.1.1 Agitate each dilution to resuspend material that may have settled during preparation. Plating should be carried out rapidly after preparing the dilutions.

7.4.1.2 Using a sterile pipette, transfer 1 mL of sample to be tested and/or 1 mL of its serial dilutions to a sterile petri dish. Samples should be plated in duplicate.

7.4.1.3 Pour 15–20 mL of agar, tempered to 44 - 47°C, into Petri dish directly on top of sample. The plate should be poured not more than 15 min after the dilution preparation.

7.4.1.4 Swirl Petri dish to homogenize sample and agar.

7.4.1.5 Allow plate to solidify.

7.4.2 Incubation

7.4.2.1 Incubate plate in the inverted position:

- 37°C for 18 - 24 h for simultaneous enumeration of *E. coli* and coliforms (= total coliforms)

AND/OR

- 44°C for 18 - 24h for enumeration of *E. coli*

7.4.2.2 It is recommended to incubate in stacks of no more than three plates. Avoid excessive crowding of plates in order to permit rapid equilibration of plates with incubator temperature.

7.5. Counting Colonies

7.5.1 After incubation, proceed with counting typical colonies on the plate.

7.5.2 Count all blue to green colonies as coliform bacteria (other than *E. coli*). If required, proceed with confirmation of presumptive positive coliforms.

7.5.3 Count all violet to pink colonies as *E. coli*. If required, proceed with confirmation of presumptive positive *E. coli*.

7.5.4 Proceed to Section 7.7 for enumeration calculations.

7.6 Confirmation (Optional)

7.6.1 For non-*E. coli* coliforms: Select two suspect non-*E. coli* coliform colonies as described in Section 7.5.2 and inoculate two tubes containing Lauryl Sulfate Tryptose (LST) broth. Incubate tubes at 35°C and examine for gas formation at 24 and 48 h.

7.6.2 Using an inoculating loop, transfer a loopful of broth from gassing tubes to Brilliant Green Lactose Bile broth (BGLB broth) and incubate at 35°C for 48 ± 2 h.

- 7.6.3 Colonies producing growth and gas in LST and BGLB are considered confirmed for non-*E. coli* coliform bacteria. The number of non-*E. coli* coliforms is determined by multiplying the percent of colonies confirmed as positive by the original RAPID'*E. coli* 2 count. The confirmed number of colonies is then multiplied by the appropriate dilution factor.
- 7.6.4 For *E. coli* coliforms: Select two suspect *E. coli* coliform colonies as described in Section 7.5.3 and inoculate two tubes containing LST broth. Incubate tubes at 35°C and examine for gas formation at 24 and 48 h.
- 7.6.5 Using an inoculating loop, transfer a loopful of broth from gassing tubes to EC (*E. coli*) broth and incubate at 45.5°C for 48 ± 2 h.
- 7.6.6 Streak a loopful of broth from gas-positive tubes onto Levine's Eosin Methylene Blue agar (L-EMB) and incubate at 35°C for 24 ± 2 h. Select a minimum of 2 typical colonies with dark purple center and green metallic sheen from each L-EMB plate for confirmation. For each colony selected, perform a Gram stain and inoculate an additional tube of LST broth- incubate at 35°C for 48 ± 2 h, and examine for gas formation. Conduct additional biochemical tests as follows: indole, methyl red, Vogues-Proskauer, and citrate tests. These biochemical tests may be performed by using a rapid ID kit that encompasses the tests specified above.
- 7.6.7 Confirmed *E. coli* isolates are Gram-negative, exhibit growth and gas in LST, are indole and methyl red positive, and Vogues-Proskauer and citrate negative.
- 7.6.8 Total number of *E. coli* colonies is determined by multiplying the percent of isolates confirmed as positive by the original RAPID'*E. coli* 2 *E. coli* count. The confirmed number of *E. coli* is then multiplied by the appropriate dilution factor. Since *E. coli* is a coliform bacteria, enumeration of total coliforms is achieved by adding the number non-*E. coli* coliforms to the number of *E. coli*.

7.7 Recording Results

- 7.7.1 To achieve enumeration of total number of coliform bacteria, add the number of blue to green and violet to pink colonies and multiply by the appropriate dilution factor.
- 7.7.2 Select only dishes containing fewer than 300 total coliform colonies in all.
- 7.7.3 Calculate the number of CFU per g or mL by multiplying the number of colonies present by the appropriate dilution factor.

Note: *E. coli* only CFU/g or mL = All violet to pink colonies x dilution factor.

Other coliforms CFU/g or mL = All blue to green colonies x dilution factor.

Total coliform CFU/g or mL = All blue to green and violet to pink colonies x dilution factor

Note: Whenever possible, the actual count should be reported. Other options, considerations or helpful hints include the following:

A: If no colonies are recovered in the 1 in 10 dilution, report as < 10 CFU/g or mL.

B: If more than 300 colonies in all appear, the food sample may be retested using a higher dilution (i.e., 1 in 50 or 1 in 100) to enable a more accurate enumeration or reported as >300 x the dilution factor.

8. References

- 8.1 Bio-Rad Laboratories – RAPID *E. coli* 2 agar Technical Sheet. Available at <http://www.bio-rad.com/en-ca/product/rapide-coli-2-medium>.
- 8.2 *Official Methods of Analysis* (1995) 16th Ed., AOAC INTERNATIONAL, Gaithersburg, MD, Method 966.24

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