Health Products and Food Branch

Ottawa

Detection of *Listeria* spp. from Environmental Surfaces and Heat Processed Ready to Eat Meat  
and Poultry Using iQ-Check™ *Listeria* spp. Real-Time PCR Test Kit

Microbiological Methods Committee  
Microbiology Evaluation Division  
Bureau of Microbial Hazards, Food Directorate,  
Health Products and Food Branch, Health Canada  
Postal Locator: 2204E  
Ottawa, Ontario K1A 0K9

Contact the Microbiological Methods Committee: [mmc-cmm@hc-sc.gc.ca](mailto:mmc-cmm@hc-sc.gc.ca)

1. Application

This method is applicable to the rapid detection of *Listeria* spp. to determine compliance with the requirements of Section 4 and 7 of the *Food and Drugs Act*. This method has been validated for use in environmental samples and the heat processed food type of the ready-to-eat meat and poultry food category. This revised method replaces MFLP-39, dated April 2011.

**Note:** While this method is only approved for environmental surfaces and heat processed ready-to-eat meat and poultry, it is assumed that this method could be used with other foods. To ensure the method is fit for purpose for commodities outside the application, it is imperative that other commodities be properly validated following the criteria in the *Compendium of Analytical Methods*. It is requested that these validation data be forwarded to the Microbiological Methods Committee so the Application Section can be expanded to include these new foods if the data comply with MMC requirements (refer to Development of Methods in Volume 1 of the *Compendium of Analytical Methods*).

2. Description

The iQ-Check™ system is a real-time PCR based assay that involves simultaneous amplification and detection of the target organism after a single overnight enrichment. The *Listeria* spp. test detects all species of *Listeria*, including *L. grayi*. PCR is run after a single 24 hr enrichment in *Listeria* Special Broth (LSB). This broth was developed by Bio-Rad and specially formulated to meet the growth requirements of *Listeria*. The nutritive nature of LSB combined with the sensitivity and specificity of the iQ-Check assay eliminate the need for a secondary enrichment. The assay should only be used by trained laboratory personnel who follow standard Good Laboratory Practices.

### 3. Principle

The iQ-Check *Listeria* spp. kit is a test based on gene amplification and detection by real-time PCR. Ready-to-use PCR reagents contain DNA primers and a fluorescent DNA hybridization probe specific for *Listeria* spp., as well as DNA polymerase and nucleotides. All reagents necessary for the test are contained in the kit. PCR is a technique used to generate many copies of target DNA. During the PCR reaction, several cycles of heating and cooling allow DNA denaturation, by heat, followed by primers binding to the target region. The DNA polymerase then uses these primers and deoxynucleotide triphosphates (dNTPs) to extend the DNA, creating copies of the target DNA. These copies are called amplicons. In real-time PCR, specific oligonucleotide probes are used to detect the DNA during the amplification, by hybridizing to the amplicons. These probes are linked to a fluorophore which fluoresces only when hybridized to the target sequence. In the iQ-Check *Listeria* spp. kits, carboxyfluorescein (FAM) is the fluorophore linked to the probe hybridizing to the *Listeria* spp. specific DNA sequence. In the absence of target DNA, no fluorescence will be detected, and the sample is determined to be negative by the iQ-Check software. As the amount of amplicons increases with each round of amplification, fluorescence intensity also increases. During each PCR cycle, at the annealing step, the real-time PCR system measures this fluorescence and the associated software plots the fluorescence intensity versus number of cycles. This method allows a simple determination of the presence of *Listeria* spp. in a sample. To monitor for a successful DNA amplification in each reaction well, a synthetic DNA "internal control" is included in the reaction mix. This control is amplified with a specific probe at the same time as the *Listeria* spp. target DNA sequence, and detected by a second fluorophore.

### 4. Definition of Terms

See [Appendix A of Volume 3](#).

### 5. Collection of Samples

See [Appendix B of Volume 3](#).

### 6. Materials and Special Equipment

**Note:** The Laboratory Supervisor must ensure that the analysis described in this method is carried out in accordance with the International Standard referred to as "ISO/IEC 17025:2005 (or latest version): General Requirements for the Competence of Testing and Calibration Laboratories".

**Note:** It is the responsibility of the laboratory to ensure equivalency if any variations of the media formulations listed here are used (either product that is commercially available or made from scratch). Please forward equivalency data to the Editor of *Compendium of Analytical Methods* for consideration of modification of this method.

#### 6.1 iQ-Check *Listeria* spp. kit (Bio-Rad catalog # 357-8113)

- 6.1.1 Reagent A - Lysis reagent, 1 bottle (20 mL), contains magnetic stirbar
- 6.1.2 Reagent B - Fluorescent probes, 1 tube (0.55 mL)
- 6.1.3 Reagent C - Amplification mix, 2 tubes (2.2 mL each)
- 6.1.4 Reagent D - PCR negative control, 1 tube (0.5 mL)
- 6.1.5 Reagent E - PCR positive control, 1 tube (0.25 mL)
- 6.1.6 Reagent F - Lysis beads, 1 bottle (17.6g)

## 6.2 Supplies and reagents

- 6.2.1 Combitips pipette or equivalent repeat pipette
- 6.2.2 Combitips tips
- 6.2.3 Conical screw cap sterile tubes, 1.5 mL (for tube extraction only)
- 6.2.4 Dey-Engley (D/E) neutralizing broth for sponges and swabs
- 6.2.5 Deepwell microplates (for deepwell extraction only)
- 6.2.6 Enrichment medium - LSB (Bio-Rad catalog # 355-5703 [225mlx6], 356-4706 [500g])
- 6.2.7 Environmental sponges
- 6.2.8 Environmental swabs
- 6.2.9 Low-profile 0.2 mL PCR 8-tube strips without caps, thin-walled, white (Bio-Rad, catalog # TLS-0851XTU)
- 6.2.10 Low-profile unskirted PCR plates, 48 or 96 well, white (Bio-Rad catalog# MLL-4851XTU, MLL-9651XTU)
- 6.2.11 Micropipettes - 20 µL, 200 µL and 1000 µL
- 6.2.12 Optical Flat 8-cap strips, ultraclear (Bio-Rad catalog # TCS-0803XTU)
- 6.2.13 Pre-pierced sealing film (Bio-Rad catalog # 360-0040) (for deepwell extraction only)
- 6.2.14 Powder-free gloves
- 6.2.15 Sterile filter tips - adaptable to 20 µL, 200 µL and 1000 µL micropipettes
- 6.2.16 Stomacher bag with incorporated filter

## 6.3 Equipment

- 6.3.1 Bench top centrifuge - maximum 10,000 – 12,000 × g, for 1.5 mL tubes (for tube extraction only)
- 6.3.2 Cell disruptor (Bio-Rad catalog # 359-1456) (for tube extraction only)
- 6.3.3 Dry heat block, for 1.5 mL tubes - capable of maintaining 95 – 100°C (for tube extraction only)
- 6.3.4 Incubator - capable of maintaining temperatures of 30 ± 1°C for sample incubation
- 6.3.5 Magnetic stir plate
- 6.3.6 Real-time PCR system – CFX96™ real-time PCR detection system (Bio-Rad cat. # 360-0037) or MiniOpticon™ real-time PCR detection system (Bio-Rad catalog # 359-1592)
- 6.3.7 Stomacher, blender or equivalent
- 6.3.8 Thermomixer (for deepwell extraction only)
- 6.3.9 Vortex apparatus

**Note:** It is the responsibility of each laboratory to ensure that the temperatures of the incubators or water baths are maintained at the recommended temperatures. Where 35°C is recommended in text of the method the incubator may be at 35 ± 1.0°C. Similarly, lower temperatures of 30 or 25 may be ± 1.0°C. However, where higher temperatures are recommended, such as 43 or 45.5°C, it is imperative that the incubators or water baths be maintained within 0.5°C due to potential lethality of the higher temperatures on the microorganism(s) being isolated.

## 7. Procedure

The test shall be carried out in accordance with the following instructions:

### 7.1 Handling of sample units

**Note:** It is important that the D/E neutralizing broth listed in 6.2 is used to moisten sponges. The use of other neutralizing buffers can cause PCR inhibition.

- 7.1.1 During storage and transport, with the exception of shelf-stable products, keep the sample units refrigerated. Keep sample units of frozen products frozen. Thaw frozen samples in a refrigerator, or under time and temperature conditions which prevent microbial growth or death.
- 7.1.2 Add D/E neutralizing broth to sterile sponges and swabs to moisten.
- 7.1.3 The protocol of MFLP-41 should be followed for obtaining and transporting swabs.
- 7.1.4 Analyze the sample units as soon as possible after receipt at the laboratory.

### 7.2 Preparation for analysis

- 7.2.1 Turn on heating block or thermomixer in advance to reach 95 – 100°C.
- 7.2.2 Turn on the thermal cycler and computer.
- 7.2.3 Carefully pour all the contents from Reagent F (lysis beads) into Reagent A (lysis reagent).

### 7.3 Preparation of sample

- 7.3.1 Temper LSB to  $30 \pm 2^\circ\text{C}$  before use.
- 7.3.2 To ensure a representative analytical unit, agitate liquids or free flowing materials until the contents are homogeneous. If the sample unit is a solid, obtain the analytical unit by taking a portion from several locations within the sample unit.
- 7.3.3 Add 25 g or mL of the food (analytical unit) to 225 ml of LSB. For composite samples, analytical units may be combined up to 125g or mL (e.g., 125 g or mL in 1125 mL of LSB). If alternate analytical units are required, maintain a ratio of 1 part sample material to 9 parts LSB.
- 7.3.4 Add the environmental sponge or large swabs to 100 mL of LSB or composite up to 10 sponges with 100 mL of LSB **per sponge**. Place smaller environmental swabs (e.g., cotton tip) in 10 mL of LSB or composite up to 10 swabs with 10 mL LSB **per swab**.
- 7.3.5 Blend, stomach or vortex as required for thorough mixing.
- 7.3.6 Incubate food sample or environmental sponges/swabs in LSB at  $30 \pm 1^\circ\text{C}$  for  $25 \pm 1$  h.

### 7.4 DNA extraction – tube format

- 7.4.1 Collect 100  $\mu\text{L}$  of enriched sample into a 1.5 mL conical tube with screw cap. Avoid transferring food particles. Store the remaining enrichment broth at  $2 - 8^\circ\text{C}$  for later use in the confirmation of any iQ Check™ *Listeria* spp. positive assay results.
- 7.4.2 Place final lysis reagent (reagent + beads) on a magnetic stir plate and stir on medium speed for a minimum of 1 min.

- 7.4.3 Add 100 µL of the final lysis reagent to the sample tube. Mix the solution by pipetting up and down in the tube a minimum of 3 times.
- 7.4.4 Place tube in the cell disruptor for 3 ± 1 min.
- 7.4.5 Place tube in the heat block at 100°C for 15 min.
- 7.4.6 Vortex at high speed.
- 7.4.7 Centrifuge at 10,000 – 12,000 × g for 2 min.

## 7.5 DNA extraction – deepwell format

- 7.5.1 Add 100 µl of the lysis reagent to each well of a deepwell microplate.
- 7.5.2 After sample enrichment, collect a 100 µl aliquot of each sample and add into the deepwell microplate. Avoid transferring food particles..
- 7.5.3 Seal with pre-pierced sealing film.
- 7.5.4 Place in the thermomixer instrument set at 95 – 99°C at 1,300 rpm for 15 – 20 min.
- 7.5.5 Remove plate from thermomixer and allow to cool in the refrigerator for 10 minutes.

## 7.6 PCR

- 7.6.1 Prepare PCR mix combining the amplification solution (Reagent C) and the fluorescent probes (Reagent B) according to the chart in the package insert.
- 7.6.2 Add 45 µL of PCR mix to each well of a 48 or 96 well PCR plate, pipetting carefully to avoid bubbles.
- 7.6.3 Add 5 µL of each sample, or control, to the appropriate well.
- 7.6.4 One positive control (Reagent E) and one negative control (Reagent D) are run with each assay. Controls are provided as ready-to-use solutions in the kit.
- 7.6.5 Seal the plate with the optical caps.
- 7.6.6 Place the plate into the reaction module of the thermal cycler and close the lid.
- 7.6.7 Verify the plate setup and protocol selection are correct.
- 7.6.8 Click on the “run” button in the software program.
- 7.6.9 A window will appear allowing you to save your data in a specific location. Once “save” is selected, the PCR protocol is started.
- 7.6.10 A report will be generated at the end of the run assigning each sample a positive or negative result.

## 7.7 Confirmation of Presumptive Positive Results

Using the refrigerated enrichment broth, proceed with plating and confirmation steps of a cultural method (*i.e.*, MFHPB-07 or MFHPB-30), as appropriate. Use at least one chromogenic agar listed in the cultural method in addition to Oxford agar.

**Note:** For the heat processed food type of the ready-to-eat meat and poultry food category, if the results indicate that *Listeria* spp. have not been detected, it is understood that the presence of *L. monocytogenes* has also not been detected, since no species belonging to the *Listeria* genus were detected. However, if *Listeria* spp. are detected in heat processed ready-to-eat meat and poultry, regardless of the amount of colony screening or additional confirmation steps done, these results cannot be used to state that *L. monocytogenes* is absent. For heat processed ready-to-eat meat and poultry, report all species of *Listeria* that are identified.

## 8. References

- 8.1 Bio-Rad. iQ-Check™ *Listeria* spp. Kit User Guide 357-8113. Test for the real-time PCR detection of *Listeria* spp. in food and environmental samples.
- 8.2 Microbiological Methods Committee. 2010. [MFLP-41: Environmental sampling for the detection of microorganisms.](#) Volume 3. *The Compendium of Analytical Methods*.
- 8.3 Pagotto, F., Hébert, K., J. Farber. 2011. [MFHPB-30: Isolation of \*Listeria monocytogenes\* and other \*Listeria\* spp. from foods and environmental samples.](#) Volume 2. *The Compendium of Analytical Methods*.
- 8.4 Warburton, D., Boville, A., Pagotto, F., Daley, E. and Chow, C. 2011. [MFHPB-07: The Isolation of \*Listeria monocytogenes\* and other \*Listeria\* spp. from foods and environmental samples using palcam broth.](#) Volume 2. *The Compendium of Analytical Methods*.

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