



Update to EN 61010-1 Amendment Statement of Operating Conditions

Catalog #	Description
1652660	Gene Pulser Xcell Total System
1652661	Gene Pulser Xcell Eukaryotic System
1652662	Gene Pulser Xcell Microbial System
1652666	Gene Pulser Xcell Main Unit, 100/240 V
1652667	Gene Pulser Xcell CE Module, 25–3,275 μ F
1652668	Gene Pulser Xcell PC Module, 50–1,000 ohm

For research purposes only.

Specifications

See the following user manual updates according to the guidelines in IEC/EN/UL 61010-1/A1, section 5.4.2:

Environmental Conditions

Operating temperature	0–35°C, indoor use
Humidity	0–95%, in the absence of condensation
Altitude	Up to 2,000 m
Pollution degree	2

Input Power

Supply voltage fluctuations	\pm 10% (unless otherwise specified), Maximum current 2 A Maximum power 240 W
Overvoltage category	II

BIO-RAD and GENE PULSER are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner.